ASCE - Telescopic Elbow ROM

EN ELBOW ORTHOSIS

DE ELLENBOGENORTHESE

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IT ORTESI PER GOMITO

NL ELLEBOOGORTHESE

DK ALBUEORTOSE







FI KYYNÄRNIVELORTOOSI

SE ARMBÅGSORTOS

NO ORTOSE, ALBUE

PL ORTEZA ŁOKCIA

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ENGLISH

Product Description

A post op elbow brace to provide controlled range of motion at the elbow joint. Telescoping bars can be adjusted to customise brace length.

Usage Information (See diagram)

- Unfasten the straps and apply the elbow brace in position as shown, ensuring it is aligned centrally with the ROM hinge at the elbow. Note: it may be required to shorten the straps - simply remove the Y-Tab, trim strap to correct length and reattach the Y-Tab in position.
- If required, press the button to extend the brace length to accommodate the patient's arm.
- The supplied sling strap can be added if required. With the injured arm supported, loop the wrist strap through the D-Ring as shown. Position the strap over the shoulder and under the opposite arm as shown, then attach the clip onto the D-Ring.

Hinge Instructions:

- Extension limit setting can be selected between -10° (Hyperextension) and 70°. Simply pull and slide the blue stopper to the desired position.
- Flexion limit settings can be selected between -10° and 120°. Simply pull and slide the other blue stopper to the desired position.
- The hinge can be locked by sliding the lock button to the locked position at any of the 5 settings: -10° (Hyperextension), 0° (neutral), 10°, 20°, 30° of flexion.

Indications for Use

- Following surgery of the elbow joint.
- For immobilisation after elbow joint injuries and subsequent mobilisation.
- Elbow ligament injuries.

Contraindications

- Do not use over open wounds.

Warnings and Precautions

We recommend the initial fitting of this brace be conducted by a suitably qualified healthcare professional who will advise the period of use.

- Carefully read all instructions and warnings prior to use.
- Follow all instructions to ensure proper performance of the brace.
- Do not use if liniments, ointments, gels, creams or any other substances have been applied to the affected area.
- Do not re-use for another patient, doing so risks cross infection and can compromise product integrity.
- Perform regular skin and circulation checks, especially patients with diabetes, vascular deficiencies and neurological conditions.
- Should any adverse reactions occur, please cease use and contact your healthcare professional or provider of this product.
- The durability of the brace may be compromised by certain factors, e.g., objects with sharp edges or damage to the hook and loop fasteners. To prevent this the hook and loop fasteners should always be fastened when the brace is not being worn or when being washed.

Washing and Care Instructions

- Close over all straps.
- Hand wash sling in cool water using mild detergent.
- Rinse and dry flat.
- Do not spin or tumble dry.
- Do not use fabric conditioner at any time.
- Check for any degradation that may affect use.
- Periodically check fasteners are operating securely.

Product Family Composition

Nylon, Aluminium, Natural Rubber, Stainless Steel, Steel, Polyester, Cotton, PVC, EVA, Glass Fibre, Silicone, POM.

Storage and Transport Conditions

Store in a cool, dry place out of direct sunlight, in the original packaging.

Recycling and Disposal

Packaging and constituent parts should be recycled or safely disposed in accordance with local or national laws.

Serious Incident

Report any serious incident to the manufacturer and the competent authority of the EU Member state, or country in which you reside.

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