

# Dorso Rigid 35 (Hewittbrace)



#### Ergonomic design

Well-considered design combines wearer comfort and effective performance.



#### Environment

We take into account the patient as well as the environment.



#### Comfort



The quality of our products shows itself in the level of comfort experienced.



The Dorso Rigid 35 (Hewitt Brace) has the following product features:

- Optimized fit and control
- Lightweight telescopic frame made of aluminum alloy
- Back cushions adjustable horizontally or vertically
- Adjustable lateral pads
- Hinged pelvic band and coccyx pad
- Spring-assisted pelvic band hinges or locks at 10° or 18°
- Coccyx pad self-adjusts for comfortable control
- User-friendly locking for optimal support
  - Effortless to put on and tighten for improved compliance
  - Safe and secure closure
  - Quick-connect strap button with hand loop for easy application

Available in the following sizes:

Article number	Article description			Size
089001/S	Dorso Rigid 35 (Hewittbrace)	60 - 75 cm	34 - 32 cm	S
089001/MS	Dorso Rigid 35 (Hewittbrace)	75 - 90 cm	34 - 41 cm	M-S
089001/M	Dorso Rigid 35 (Hewittbrace)	75 - 90 cm	36 - 46 cm	M
089001/LS	Dorso Rigid 35 (Hewittbrace)	90 - 105 cm	36 - 46 cm	L-S
089001/L	Dorso Rigid 35 (Hewittbrace)	90 - 105 cm	40 - 52 cm	L
089001/XL	Dorso Rigid 35 (Hewittbrace)	105 - 115	42 - 54 cm	XL