

Elbow Brace



Ergonomic design

Well-considered design combines wearer comfort and effective performance.



Environment

We take into account the patient as well as the environment.



Comfort

The quality of our products shows itself in the level of comfort experienced.



The Elbow ROM Brace is easy to adjust, providing excellent comfort and support for recovery from surgery, injury, or instability.

The range of motion can be easily adjusted using the spring-loaded quick clamps with a simple 'pull and slide' movement. The design prevents accidental changes to the flexion and extension settings.

Available in the following sizes:

Article number	Article description	Size
139503	Elbow Brace right	Universal
139505	Elbow Brace left	Universal