

# Ligaflex Rhizo Thumb Brace



#### Ergonomic design

Well-considered design combines wearer comfort and effective performance.



#### Environment

We take into account the patient as well as the environment.



#### Comfort

The quality of our products shows itself in the level of comfort experienced.



This thumb brace for immobilization of the thumb in an analgesic position perfectly adapts to the morphology thanks to a metal wire frame with a spring effect. The Thermytex material on the inside is breathable and antibacterial. Two Velcro closures simplify the application of the product. Immobilization of the thumb and the base of the hand while preserving hand function.

#### Indications:

- Trauma to the base of the hand or thumb (CMC and/or MCP-joints)
- Rehabilitation after surgery of the base of the hand or thumb
- CMC/ MCP osteoarthritis
- Rheumatoid arthritis in the hand
- Immobilization of the trapeziometacarpal joint in cases of early-stage rhizarthrosis
- Pain in the thumb base joint
- Osteoarthritis around the scaphoid bone
- Postoperative care
- Application of the prosthesis
- Trapeziectomy
- Surgery of metacarpal fractures at the base of MC I: Rolando or Bennett fractures
- Mild thumb sprains

#### Available in the following sizes:

Right	Left	Article description	Wrist circumference	Size
GMI083	GMI093	Ligaflex Rhizo Thumb Brace	13 - 17 cm	1
GMI085	GMI095	Ligaflex Rhizo Thumb Brace	17 - 23 cm	2