

Ligaflex Classic open

Wrist Brace



Ergonomic design

Well-considered design combines wearer comfort and effective performance.



Environment

We take into account the patient as well as the environment.



Comfort

The quality of our products shows itself in the level of comfort experienced.



The Classic Open Wrist Brace is an economical, ready-to-use brace for immobilizing the wrist. The wrist brace features a removable aluminum palmar splint and three Velcro straps for an ideal fit. Due to its smart design, this brace can be applied to both the left and right wrists. Additionally, this brace is easily adjustable, requiring only two sizes.

Indications:

- Post-traumatic wrist pain
- Post-operative wrist pain
- Carpal tunnel syndrome (CTS)

Available in the following sizes:

Item Name	Article description	Wrist circumference	Size
243702-1	Ligaflex Classic open Wrist Brace	13 - 17 cm	1
243702-2	Ligaflex Classic open Wrist Brace	17 - 23 cm	2