

Lumbosacral Brace Low

(27 cm)



Ergonomic design

Well-considered design combines wearer comfort and effective performance.



Environment

We take into account the patient as well as the environment.



Comfort

The quality of our products shows itself in the level of comfort experienced.



The brace helps to recover as quickly and effectively as possible. The Lumbosacral Brace Low has a height of 27 cm and is anatomically designed to provide support to the lumbosacral area. The brace features a wide, easy-to-fasten Velcro closure at the front. The four adjustable Velcro straps allow for better regulation and adjustment of tension in the back area.

Comfortable – Unisex – Hygienic: Soft cotton in contact with the skin, made from breathable woven material.

Indications:

- Lumbago
- Sciatica
- Lumbar protrusions
- Lumbar osteoarthritis
- Lumbar muscle contractures
- Discopathy
- Lumbar hyper/hypo lordosis

Available in the following sizes:

Article number	Article description	Circumference of the waist in cm	Size
GM8582-S	Lumbosacral Brace Low (27 cm)	80 - 90 cm	S
GM8582-M	Lumbosacral Brace Low (27 cm)	90 - 100 cm	M
GM8582-L	Lumbosacral Brace Low (27 cm)	100 - 110 cm	L
GM8582-XL	Lumbosacral Brace Low (27 cm)	110 - 120 cm	XL
GM8582-XXL	Lumbosacral Brace Low (27 cm)	120 - 130 cm	XXL
GM8582-XXXL	Lumbosacral Brace Low (27 cm)	> 130 cm	XXXL