

Ortel C4 Vario (Miami)



Ergonomic design

Well-considered design combines wearer comfort and effective performance.



Environment

We take into account the patient as well as the environment.



Comfort

The quality of our products shows itself in the level of comfort experienced.



The Ortel C4Vario is a universal cervical collar. It is fully adjustable for optimal support. Extra support for the cervical spine is provided by two rigid plates. The occipital, chin, high back, and sternal supports ensure complete immobilization. The brace is adjustable in height from 9 to 16 cm. The front opening allows access to the trachea. The brace is comfortable due to its removable and washable lining and is radiolucent.

Indications:

- Severe neck pain
- Stable fractures
- Osteoarthritis

Available in the following sizes:

Article number	Article description	Color	Size
T492901	Ortel C4Vario (Miami)	Gray	Universal