

Ortho E-Back Brace

(H33 cm)



Ergonomic design

Well-considered design combines wearer comfort and effective performance.



Environment

We take into account the patient as well as the environment.



Comfort

The quality of our products shows itself in the level of comfort experienced.



The Ortho E-Back Brace features an anatomical design that ensures a perfect fit for the patient. The semi-rigid textile in the paravertebral area provides maximum stability for the spine, while the additional front straps effectively enhance the supportive effect. The number of straps varies depending on the height of the corset (to adjust tension in different parts of the back). The stays on the back side can be removed to adapt the support to the patient's needs and different rehabilitation phases. The innovative bacteriostatic Radyarn fabric guarantees maximum hygiene, comfort, breathability, and prevents unpleasant odors.

Indications:

- Acute or chronic lower back pain and/or lumbar sciatica
- Moderate spondyloarthritis
- Spinal muscle contractures
- Lumbosacral tendinopathies
- Disc herniation
- Lumbosacral kyphosis
- Mild osteolysis
- Postoperative after a hernia surgery, posterior decompression

Available in the following sizes:

Article number	Article description	Circumference of the waist in cm	Size
GM8553.33-S	Ortho E-Back Brace (H33 cm)	80 - 90 cm	S
GM8553.33-M	Ortho E-Back Brace (H33 cm)	90 - 100 cm	M
GM8553.33-L	Ortho E-Back Brace (H33 cm)	100 - 110 cm	L
GM8553.33-XL	Ortho E-Back Brace (H33 cm)	110 - 120 cm	XL
GM8553.33-XXL	Ortho E-Back Brace (H33 cm)	120 - 130 cm	XXL
GM8553.33-XXXL	Ortho E-Back Brace (H33 cm)	> 130 cm	XXXL