

P-FIX + Wrist Splint



Ergonomic design

Well-considered design combines wearer comfort and effective performance.



Environment

We take into account the patient as well as the environment.



Comfort

The quality of our products shows itself in the level of comfort experienced.



The P-FIX+ brace features four anatomically shaped, moldable, and removable splints. It is equipped with three velcro straps to adjust the closure according to the size of the patient's forearm. The brace is fully openable, making it easy to apply to painful or recently operated wrists. The bacteriostatic Radyarn® textile ensures improved hygiene, comfort, and breathability.

Indications:

- **Traumatic conditions:** sprains, carpal and metacarpal fractures, tendonitis, De Quervain's tenosynovitis, thumb traumas.
- **Degenerative conditions:** osteoarthritis, rhizartrosis.
- **Rehabilitation phase:** post-traumatic and postoperative relief.

Available in the following sizes:

Right	Left	Article description	Size
GM8134-RS	GM8134-LS	P-FIX+ Thumb and Wrist Splint	S
GM8134-RM	GM8134-LM	P-FIX+ Thumb and Wrist Splint	M
GM8134-RL	GM8134-LL	P-FIX+ Thumb and Wrist Splint	L