

Wrist Splint



Ergonomic design

Well-considered design combines wearer comfort and effective performance.



Environment

We take into account the patient as well as the environment.



Comfort

The quality of our products shows itself in the level of comfort experienced.



The Wrist Splint is a ready-made brace for immobilization of the wrist. The Wrist Splint features a removable aluminum palmar stay and three Velcro straps for an ideal fit. The soft thumb opening ensures a high level of comfort.

Available in the following sizes:

Right	Left	Width of the hand at the M.C.P. joint	Size
GMI011	GMI021	5,0 - 6,3 cm	XX-Small*
GMI012	GMI022	6,3 - 7,5 cm	X-Small
GMI013	GMI023	7,5 - 8,8 cm.	Small
GMI014	GMI024	8,8 - 10,0 cm.	Medium
GMI015	GMI025	10,0 - 11,3 cm.	Large
GMI016	GMI026	11,3 - 12,5 cm.	X-Large

*Size XX-Small is used for (slightly bigger) children, approximately ages between 9-12 years.