

Moldable Wrist Brace



Ergonomic design

Well-considered design combines wearer comfort and effective performance.



Environment

We take into account the patient as well as the environment.



Comfort

The quality of our products shows itself in the level of comfort experienced.



The Moldable Wrist Brace is 25 cm long and is used with an airtight packaged polyester splint. This splint can be applied in any desired position after opening the packaging and immersing it in water. The polyester splint is provided separately and can be inserted into the wrist brace using a specially designed tweezer.

Available in the following sizes:

Right	Left	Article description	Wrist circumference	Size
GMI053-I	GMI063-I	Moldable Wrist Brace	7,5 - 8,8 cm	S
GMI054-I	GMI064-I	Moldable Wrist Brace	8,8 - 10,0 cm	M
GMI055-I	GMI065-I	Moldable Wrist Brace	10,0 - 11,3 cm	L
GMI056-I	GMI066-I	Moldable Wrist Brace	11,3 - 12,5 cm	XL

