

Moldable Wrist Brace with Thumb



Ergonomic design

Well-considered design combines wearer comfort and effective performance.



Environment

We take into account the patient as well as the environment.



Comfort

The quality of our products shows itself in the level of comfort experienced.



The Moldable Wrist Brace with Thumb consists of a wrist brace (25 cm long) and an airtight packaged polyester splint that, after opening and immersion in water, serves as a stay and can be positioned in any desired position. The polyester splint is inserted into the wrist brace using a specially designed tweezers.

Available in the following sizes:

Right	Left	Article description	Width of the hand at the level of the M.C.P. joint	Size
GMI033-I	GMI043-I	Moldable Wrist Brace with Thumb	7,5 - 8,8 cm	S
GMI034-I	GMI044-I	Moldable Wrist Brace with Thumb	8,8 - 10,0 cm	M
GMI035-I	GMI045-I	Moldable Wrist Brace with Thumb	10,0 - 11,3 cm	L
GMI036-I	GMI046-I	Moldable Wrist Brace with Thumb	11,3 - 12,5 cm	XL

Tweezers moldable wrist brace

GMI070	30 cm
--------	-------