

# Stabilizing Shoulder Brace



## Ergonomic design

Well-considered design combines wearer comfort and effective performance.



## Environment

We take into account the patient as well as the environment.



## Comfort

The quality of our products shows itself in the level of comfort experienced.



This stabilizing shoulder brace helps ensure a fast and effective recovery. The straps fasten upwards to keep the humerus in its natural position (reducing pain and recovery time). It features a silicone inner arm ring for maximum grip and better arm support.

An elastic extra strap pushes the shoulder backward in cases of anterior shoulder dislocations. The brace is comfortable: lightweight and breathable. Hygienic: 100% cotton in contact with the skin.

## Indications:

- Shoulder dislocations and chronic shoulder subluxation
- Prevention of dislocations
- Osteoarthritis of the shoulder
- Consequences of previous dislocations of the scapulohumeral joint
- Post-operative care
- Stabilization of the glenohumeral joint
- Conservative treatment for neurological shoulder issues

## Available in the following sizes:

| Article number | Article description        | Circumference of biceps | Size |
|----------------|----------------------------|-------------------------|------|
| GM8407-S       | Stabilizing Shoulder Brace | 24-28 cm                | S    |
| GM8407-M       | Stabilizing Shoulder Brace | 28-32 cm                | M    |
| GM8407-L       | Stabilizing Shoulder Brace | 32-36 cm                | L    |
| GM8407-XL      | Stabilizing Shoulder Brace | 36-40 cm                | XL   |
| GM8407-XXL     | Stabilizing Shoulder Brace | 40-44 cm                | XXL  |