

Elastic Splint Straps



Ergonomic design

Well-considered design combines wearer comfort and effective performance.



Environment

We take into account the patient as well as the environment.



Comfort

The quality of our products shows itself in the level of comfort experienced.



The Splint Straps (standard white) are made from woven polyamide. A portion of the loop strap is connected to a portion of the hook strap using ultrasonic welding. By adhering the hook section to the desired material and then pressing the loop strap against it, you achieve a reliable and secure closure that can be adjusted each time.

The Splint Straps are available in various sizes. Colors are available upon request.

Available in the following sizes:

Article number	Article description	Dimensions in cm	Color
GM941	Elastic Splint Straps	2,5 x 30	White
GM942	Elastic Splint Straps	2,5 x 40	White
GM943	Elastic Splint Straps	2,5 x 50	White
GM944	Elastic Splint Straps	2,5 x 60	White
GM945	Elastic Splint Straps	2,5 x 70	White
GM946	Elastic Splint Straps	2,5 x 80	White
GM948	Elastic Splint Straps	2,5 x 100	White
GM953	Elastic Splint Straps	5,0 x 50	White
GM954	Elastic Splint Straps	5,0 x 60	White
GM958	Elastic Splint Straps	5,0 x 100	White