

Genu Dynastab Knee Brace



Ergonomic design

Well-considered design combines wearer comfort and effective performance.



Environment

We take into account the patient as well as the environment.



Comfort

The quality of our products shows itself in the level of comfort experienced.



The Genu Dynastab Knee Brace is a supportive stabilizing brace for the mildly unstable knee. The rigid stays with biaxial hinges combined with adjustable Velcro straps ensure good stability. The brace has an excellent fit, thanks to the anatomically woven design using multiple weaving techniques. The brace is available in 6 sizes and can be worn on either the left or right knee.

Indications:

- Mild to moderate partial ligament tears in the knee
- Supportive in rehabilitation for daily or sports activities
- General mild instability of the knee

Available in the following sizes:

Article number	Article description	Knee circumference	Size
237005-T1	Genu Dynastab Knee brace	32 - 34 cm	1
237005-T2	Genu Dynastab Knee brace	35 - 37 cm	2
237005-T3	Genu Dynastab Knee brace	38 - 41 cm	3
237005-T4	Genu Dynastab Knee brace	42 - 45 cm	4
237005-T5	Genu Dynastab Knee brace	46 - 49 cm	5
237005-T6	Genu Dynastab Knee brace	50 - 53 cm	6