

Junior knee immobiliser



Ergonomic design

Well-considered design combines wearer comfort and effective performance.



Environment

We take into account the patient as well as the environment.



Comfort

The quality of our products shows itself in the level of comfort experienced.



The Junior Knee Immobilizer is a comfortable splint for the immobilization of the knee. The Tri-Panel design allows for the circumference of the thigh to be individually adjusted. Various sizes are available depending on the length of the leg.

Indications:

- Immobilization brace with 4 moldable and X-ray transparent splints;
- Anti-slip system and Velcro closure to ensure the orthosis stays perfectly in place;
- Suitable for all calf and thigh circumferences thanks to 3 movable parts.

Available in the following sizes:

Article number	Inseam length	Size
I20321	50 cm - 60 cm	Size 1
I20322	60 cm - 70 cm	Size 2
I20323	70 cm - 80 cm	Size 3
I20324	80 cm - 120 cm	Size 4