

Knee Brace For Light MCL-LCL Instability



Ergonomic design

Well-considered design combines wearer comfort and effective performance.



Environment

We take into account the patient as well as the environment.



Comfort

The quality of our products shows itself in the level of comfort experienced.



This Knee Brace for Light MCL-LCL Instability provides support to the knee and allows the user to apply the desired pressure above or below the knee. The brace features two (removable) polycentric stays that provide lateral stability and support to the medial and lateral collateral ligaments. The soft, elastic, and breathable material offers high wearing comfort.

Indications:

- Traumatic pathologies: mild sprains, mild joint laxity, bursitis, joint effusion.
- During the rehabilitation phase: passive support during rehabilitation, mild protection during sports activities.

Available in the following sizes:

Article Name	Article description	Circumference in cm thigh 10 cm above the patella	Size
GM8075-S	Knee Brace For Light MCL-LCL Instability	32– 36 cm	Small
GM8075-M	Knee Brace For Light MCL-LCL Instability	36– 40 cm	Medium
GM8075-L	Knee Brace For Light MCL-LCL Instability	40– 45 cm	Large
GM8075-XL	Knee Brace For Light MCL-LCL Instability	44– 49 cm	X-Large
GM8075-XXL	Knee Brace For Light MCL-LCL Instability	48– 54 cm	XX-Large
GM8075-XXXL	Knee Brace For Light MCL-LCL Instability	54– 60 cm	XXX-Large
GM8075-XXXXL	Knee Brace For Light MCL-LCL Instability	60– 66 cm	XXXX-Large