

# Ligaflex Clavicula



#### Ergonomic design

Well-considered design combines wearer comfort and effective performance.



#### Environment

We take into account the patient as well as the environment.



#### Comfort

The quality of our products shows itself in the level of comfort experienced.



The back strap for stabilisation of the clavicle is intended for post-traumatic use. It has two shoulder straps attached to an anchoring pad on the high back. Simultaneous tightening of the straps helps stabilise the clavicle. The straps are adjustable and adaptable to the patient's shoulder width, and the product is for non-surgical clavicle fractures. For extra comfort, the band is made of a soft and lightweight fabric with padding.

#### Available in the following sizes:

Article number	Article description	Color	Chest size	Size
245003-0	Ligaflex Clavicula Junior	black	20 - 34 cm	0
245003-1	Ligaflex Clavicula	black	34 - 100 cm	I

