

Ligastep AIR walker



Ergonomic design

Well-considered design combines wearer comfort and effective performance.



Environment

We take into account the patient as well as the environment.



Comfort

The quality of our products shows itself in the level of comfort experienced.



Walker with rigid side supports to immobilize and stabilize the lower leg, ankle, and foot. The walker features an inflatable lining with an integrated pump and valve. Designed to ensure minimal dorsiflexion of the metatarsals. The non-slip roll-off sole provides patients with a natural, comfortable walking motion. The product can be worn on either foot.

Indications:

- Ankle/foot postoperative/post-traumatic immobilization
- Severe ankle/foot sprain
- Stable fracture of the lower third of the leg (long version)
- Conservative treatment of ankle/foot fractures
- Severe ankle/foot tendinopathy
- Conservative treatment of Achilles tendon rupture

Available in the following sizes:

Article number	Article description	Size	Shoe size
8403201	Ligastep AIR walker long	XS	32 - 35
8403202	Ligastep AIR walker long	S	35 - 38
8403203	Ligastep AIR walker long	M	39 - 43
8403204	Ligastep AIR walker long	L	44 - 47
8403205	Ligastep AIR walker long	XL	47 - 50

Ligastep AIR walker



Ergonomic design

Well-considered design combines wearer comfort and effective performance.



Environment

We take into account the patient as well as the environment.



Comfort

The quality of our products shows itself in the level of comfort experienced.



Walker with rigid side supports to immobilize and stabilize the lower leg, ankle, and foot. The walker features an inflatable lining with an integrated pump and valve. Designed to ensure minimal dorsiflexion of the metatarsals. The non-slip roll-off sole provides patients with a natural, comfortable walking motion. The product can be worn on either foot.

Indications:

- Ankle/foot postoperative/post-traumatic immobilization
- Severe ankle/foot sprain
- Stable fracture of the lower third of the leg (long version)
- Conservative treatment of ankle/foot fractures
- Severe ankle/foot tendinopathy
- Conservative treatment of Achilles tendon rupture

Available in the following sizes:

Article number	Article description	Size	Shoe size
8503201	Ligastep AIR walker short	XS	32 - 35
8503202	Ligastep AIR walker short	S	35 - 38
8503203	Ligastep AIR walker short	M	39 - 43
8503204	Ligastep AIR walker short	L	44 - 47
8503205	Ligastep AIR walker short	XL	47 - 50