

# Ligastep FIX walker



#### Ergonomic design

Well-considered design combines wearer comfort and effective performance.



#### Environment

We take into account the patient as well as the environment.



#### Comfort

The quality of our products shows itself in the level of comfort experienced.



This walker features rigid side supports to immobilize and stabilize the lower leg, ankle, and foot. The walker has 90° rigid uprights on both sides and is designed to ensure minimal dorsiflexion of the metatarsals. The non-slip roll-off sole provides patients with a natural and comfortable walking motion. The product can be worn on either foot and is available in both long and short versions.

#### Indications:

- Postoperative and post-traumatic immobilization of the ankle and foot
- Severe ankle and foot sprain
- Stable fracture of the lower third of the leg (long version)
- Conservative treatment of ankle and foot fractures
- Severe tendinopathy of the ankle and foot
- Conservative treatment of Achilles tendon rupture

#### Available in the following sizes:

Article number	Article description	Size	Shoe size
8003201	Ligastep FIX walker long	XS	32 - 35
8003202	Ligastep FIX walker long	S	35 - 38
8003203	Ligastep FIX walker long	M	39 - 43
8003204	Ligastep FIX walker long	L	44 - 47
8003205	Ligastep FIX walker long	XL	47 - 50

# Ligastep FIX walker



## Ergonomic design

Well-considered design combines wearer comfort and effective performance.



## Environment

We take into account the patient as well as the environment.



## Comfort

The quality of our products shows itself in the level of comfort experienced.



This walker features rigid side supports to immobilize and stabilize the lower leg, ankle, and foot. The walker has 90° rigid uprights on both sides and is designed to ensure minimal dorsiflexion of the metatarsals. The non-slip roll-off sole provides patients with a natural and comfortable walking motion. The product can be worn on either foot and is available in both long and short versions.

## Indications:

- Postoperative and post-traumatic immobilization of the ankle and foot
- Severe ankle and foot sprain
- Stable fracture of the lower third of the leg (long version)
- Conservative treatment of ankle and foot fractures
- Severe tendinopathy of the ankle and foot
- Conservative treatment of Achilles tendon rupture

## Available in the following sizes:

Article number	Article description	Size	Shoe size
8103201	Ligastep FIX walker Short	XS	32 - 35
8103202	Ligastep FIX walker Short	S	35 - 38
8103203	Ligastep FIX walker Short	M	39 - 43
8103204	Ligastep FIX walker Short	L	44 - 47
8103205	Ligastep FIX walker Short	XL	47 - 50