Ligastep PEDI Walker



Ergonomic design

Well-considered design combines wearer comfort and effective performance.



Environment

We take into account the patient as well as the environment.



Comfort

The quality of our products shows itself in the level of comfort experienced.



The Ligastep PEDI Walker is a walker for children available in three sizes. This walker features rigid side supports to immobilize and stabilize the lower leg, ankle, and foot at 90°. The walker has a rounded sole that guarantees a natural and comfortable walking pattern.

Indications:

- Postoperative and post-traumatic immobilization of the ankle and foot
- Severe ankle and foot sprain
- Stable fracture of the distal third of the leg (long version)
- Conservative treatment of ankle and foot fractures

Available in the following sizes:

Article number	Article description	Shoe size	Size
GM2651	Ligastep PEDI Walker	21 - 23	S
GM2652	Ligastep PEDI Walker	24 - 27	М
GM2653	Ligastep PEDI Walker	27 - 29	L