

# Plaster Products

Nemoa™



#### Breathable

Provides more comfort and prevents skin issues for the patient.



#### Well moldable

Moldable for optimal adjustment to anatomical differences.



#### Comfortable

The quality of a product is also reflected in the comfort for the patient.



Thin and much stronger and lighter than traditional POP. The excellent maneuverability of the bandages and optimal lamination ensure a smooth final result with rigid or semi-rigid support for the patient.

## NEMOA™ ELASTICAST RIGID POLYESTER PLASTIC CAST BANDAGE

Nemoa ElastiCast is a lightweight, strong, and slightly flexible elastic rigid plastic cast bandage. It is made from 100% polyester impregnated with polyurethane resin. Each roll is supplied in airtight packaging.

Available in the following sizes:

Article number	Article description	Size
GM521	Polyester Cast	5,0 x 360 cm
GM522	Polyester Cast	7,5 x 360 cm
GM522	Polyester Cast	10,0 x 360 cm
GM523	Polyester Cast	12,5 x 360 cm

Color White (W)

# Plaster Products

Nemoa™



#### Breathable

Provides more comfort and prevents skin issues for the patient.



#### Well moldable

Moldable for optimal adjustment to anatomical differences.



#### Comfortable

The quality of a product is also reflected in the comfort for the patient.

## NEMOA™ FLEX SEMI-RIGIDE POLYESTER KUNSTSTOFFGIPSBAND

Nemoa™ Flex is an innovative, lightweight, and strong semi-rigid plastic cast bandage. It is made from 100% polyester impregnated with polyurethane resin. The final product of Nemoa™ Flex is semi-rigid and exceptionally smooth. Each roll is supplied in airtight packaging. Nemoa™ Flex is available not only in white but also in blue, pink, and red.

Available in the following sizes:

Article number	Article description	Size
GM720	Nemoa Flex Cast	2,5 x 180 cm
GM721	Nemoa Flex Cast	5,0 x 360 cm
GM722	Nemoa Flex Cast	7,5 x 360 cm
GM723	Nemoa Flex Cast	10,0 x 360 cm
GM723	Nemoa Flex Cast	12,5 x 360 cm

Colors White (W) • Black (A) • Blue (B) • Red (R) • Pink (P)

## NEMOA™ EVER FRESH

Nemoa™ Ever Fresh is a unique, water-resistant tubular padding material suitable for use under synthetic casts. Thanks to its water resistance, this padding allows for the creation of a cast application that enables activities such as showering and (potentially) swimming.

Article number	Article description	Size
GM341	Nemoa Everfresh	5,0 x 1000 cm
GM342	Nemoa Everfresh	7,0 x 1000 cm
GM343	Nemoa Everfresh	10,0 x 1000 cm
GM344	Nemoa Everfresh	12,5 x 1000cm

# Plaster Products

Nemoa™



#### Breathable

Provides more comfort and prevents skin issues for the patient.



#### Well moldable

Moldable for optimal adjustment to anatomical differences.



#### Comfortable

The quality of a product is also reflected in the comfort for the patient.

## NEMOA™ SPLINT

The Nemoa™ Splints are available in rolls or as pre-cuts. Both are ready-to-use splints that can be applied directly without gloves, cotton, or stockinette. They are made from a 100% nonwoven polyester material impregnated with polyurethane resin. The Nemoa™ Roll splints come in a (reclosable) airtight packaging and are supplied in a convenient box. The Nemoa™ pre-cuts are individually airtight packaged and are supplied in boxes of 5.

### Available in the following sizes:

Article number	Article description	Size
GM301	Nemoa Splint Roll splints	5,0 x 450 cm
GM302	Nemoa Splint Roll splints	7,5 x 450 cm
GM303	Nemoa Splint Roll splints	10,0 x 450 cm
GM304	Nemoa Splint Roll splints	12,5 x 450 cm
GM305	Nemoa Splint Roll splints	15,0 x 450 cm

### Pre-cuts

Article number	Article description	Size
GM351	Nemoa Splint Pre-cut	5,0 x 25 cm
GM361	Nemoa Splint Pre-cut	7,5 x 30 cm
GM363	Nemoa Splint Pre-cut	7,5 x 88 cm
GM371	Nemoa Splint Pre-cut	10,0 x 38 cm
GM373	Nemoa Splint Pre-cut	10,0 x 75 cm