

# Night Splint



## Ergonomic design

Well-considered design combines wearer comfort and effective performance.



## Environment

We take into account the patient as well as the environment.



## Comfort

The quality of our products shows itself in the level of comfort experienced.



The night splint keeps the foot in dorsiflexion and allows for loading of the soft tissues of the plantar fascia and Achilles tendon. Dorsiflexion can be easily adjusted using double Velcro closures and additional wedges (3 included). Suitable for both left and right sides.

Available in the following sizes:

Article number	Size	Shoe size	
		Men	Women
U09101-S	Small	<40	<38
U09101-M	Medium	40-43	38-41
U09101-L	Large	>43	>41