

Patellar brace



Ergonomic design

Well-considered design combines wearer comfort and effective performance.



Environment

We take into account the patient as well as the environment.



Comfort

The quality of our products shows itself in the level of comfort experienced.



The patellar brace supports the patella due to its anatomical shape and exerts localized pressure on the quadriceps tendon just below the patella. The brace is comfortable, thanks to a thin strap at the back.

Indications:

- Jumper's knee
- Inflammation of the quadriceps tendon attachment
- Localized pain at the bottom of the kneecap

Available in the following sizes:

Article number	Article description	Circumference of the knee	Size
2304012	Patellar brace	26 - 50 cm	Universal

