Silistab Achillo Brace



Ergonomic designWell-considered

Well-considered design combines wearer comfort and effective performance.



Environment

We take into account the patient as well as the environment.



Comfort

The quality of our products shows itself in the level of comfort experienced.



Thanks to its combination of an elastic tricot with zones of a particular type of tricot at the ankles and a layer of silicones at the level of the Achilles tendon, this brace provides a very high level of comfortable support. The soles can be stacked to create two different heights. The brace comes with additional soles which can be used to compensate for the raised height of one leg in the other leg.

Available in the following sizes:

Article number	Article description	Ankle circumference	Size
235502-1	Silistab Achillo Brace	19 - 21 cm	I
235502-2	Silistab Achillo Brace	21 - 23 cm	2
235502-3	Silistab Achillo Brace	23 - 25 cm	3
235502-4	Silistab Achillo Brace	25 - 27 cm	4
235502-5	Silistab Achillo Brace	27 - 30 cm	5