

Ascender Telescopic ROM Knee Braces



Ergonomic design

Well-considered design combines wearer comfort and effective performance.



Environment

We take into account the patient as well as the environment.



Comfort

The quality of our products shows itself in the level of comfort experienced.



The Ascender Telescopic ROM Knee Brace is universal and easy to adjust. The universal sizing has a length range from 44 cm to 70 cm.

The ROM brace features 5 positions with a locking mechanism (drop lock) and quick-release closures for range-of-motion adjustments with a simple 'pull and slide' movement.

The design prevents unintended changes to the flexion and extension settings.

Indications:

- For comfort and support postoperatively, in case of injuries, or knee instability.

Available in the following sizes:

Article number	Article description	Length range	Size
138003	Ascender Telescopic ROM Knee Brace	44 - 70 cm	Universal