

# Ring Straps



## Ergonomic design

Well-considered design combines wearer comfort and effective performance.



## Environment

We take into account the patient as well as the environment.



## Comfort

The quality of our products shows itself in the level of comfort experienced.



The Ring Straps (standard white) are made from woven polyamide. A large portion of the loop strap is connected to a smaller portion of the hook strap using ultrasonic welding. The end of the hook section is equipped with an “easy peel,” which means that the closure can be easily opened. By threading the hook section through the buckle and then pressing it onto the loop strap, you achieve a reliable and secure closure that is easy to adjust.

The Ring Straps are available in various sizes. Colors are available upon request.

## Available in the following sizes:

Article number	Article description	Dimensions in cm
GMI00	Ring Straps	2,5 x 20
GMI01	Ring Straps	2,5 x 30
GMI02	Ring Straps	2,5 x 40
GMI03	Ring Straps	2,5 x 50
GMI04	Ring Straps	2,5 x 60
GMI05	Ring Straps	2,5 x 70
GMI06	Ring Straps	2,5 x 80
GMI07DC	Ring Straps	2,5 x 90
GMI08DC	Ring Straps	2,5 x 100
GMI10DC	Ring Straps	2,5 x 120

# Ring Straps



## Ergonomic design

Well-considered design combines wearer comfort and effective performance.



## Environment

We take into account the patient as well as the environment.



## Comfort

The quality of our products shows itself in the level of comfort experienced.



The Ring Straps (standard white) are made from woven polyamide. A large portion of the loop strap is connected to a smaller portion of the hook strap using ultrasonic welding. The end of the hook section is equipped with an “easy peel,” which means that the closure can be easily opened. By threading the hook section through the buckle and then pressing it onto the loop strap, you achieve a reliable and secure closure that is easy to adjust.

The Ring Straps are available in various sizes. Colors are available upon request.

## Available in the following sizes:

Article number	Article description	Dimensions in cm
GMI22	Ring Straps	5,0 x 40
GMI23	Ring Straps	5,0 x 50
GMI24	Ring Straps	5,0 x 60
GMI25	Ring Straps	5,0 x 70
GMI26	Ring Straps	5,0 x 80
GMI27	Ring Straps	5,0 x 90
GMI30	Ring Straps	5,0 x 120

