

White Label Comfort

Covered on both sides



Ergonomic design

Thoughtfully designed, our products combine comfort with effective functionality.



Sustainable

Optimal consideration has been given to both people and the environment.



Comfortable

The quality of a product is also reflected in the comfort it provides for the patient.



The White Label Comfort thermoplastic material is a low-temperature thermoplastic material resulting from copolymerization based on polycaprolactone and polyurethane, covered on both sides with a foam layer of 0.6 mm thickness. This unique foam layer prevents the material from sticking together and creates an air layer between the patient's skin and the thermoplastic.

The White Label Comfort thermoplastic material hardens within minutes, and adjustments can be made easily. Each sheet is individually packaged and available in various sizes, thicknesses, and types of perforations. Additionally, several pre-cut models are available (which you will find a few pages further on).

Available in the following sizes:

Article number	Type	Color	Dimensions in cm.	Thickness in mm.
WL0800	Non-Perfo	Skin color	60 x 41	0,8
WL0810	Micro-Perfo	Skin color	60 x 41	0,8
WL1600	Non-Perfo	Skin color	60 x 43	1,6
WL1610	Micro-Perfo	Skin color	60 x 43	1,6
WL2000	Non-Perfo	Skin color	60 x 44	2
WL2010	Micro-Perfo	Skin color	60 x 44	2
WL2020	Mini-Perfo	Skin color	60 x 44	2
WL2030	Multi-Perfo	Skin color	60 x 44	2
WL2500	Non-Perfo	Skin color	60 x 45	2,5
WL2520	Mini-Perfo	Skin color	60 x 45	2,5

White Label Comfort

Covered on both sides



Ergonomic design

Thoughtfully designed, our products combine comfort with effective functionality.



Sustainable

Optimal consideration has been given to both people and the environment.



Comfortable

The quality of a product is also reflected in the comfort it provides for the patient.



Available in the following sizes:

Article number	Type	Color	Dimensions in cm.	Thickness in mm.
WL2530	Multi-Perfo	Skin color	60 x 45	2,5
WL2500/9	Non-Perfo	Skin color	60 x 90	2,5
WL2520/9	Mini-Perfo	Skin color	60 x 90	2,5
WL2530/9	Multi-Perfo	Skin color	60 x 90	2,5
WL3000	Non-Perfo	Skin color	60 x 45	3
WL3020	Mini-Perfo	Skin color	60 x 45	3
WL3030	Multi-Perfo	Skin color	60 x 45	3
WL3000/9	Non-Perfo	Skin color	60 x 90	3
WL3020/9	Mini-Perfo	Skin color	60 x 90	3
WL3030/9	Multi-Perfo	Skin color	60 x 90	3
WL4000	Non-Perfo	Skin color	60 x 45	4
WL4020	Mini-Perfo	Skin color	60 x 45	4
WL4000/9	Non-Perfo	Skin color	60 x 90	4
WL4020/9	Mini-Perfo	Skin color	60 x 90	4

White Label Comfort

One Side Coated

Coated on one side with a foam layer



Ergonomic design

Thoughtfully designed, our products combine comfort with effective functionality.



Sustainable

Optimal consideration has been given to both people and the environment.



Comfortable

The quality of a product is also reflected in the comfort it provides for the patient.



The White Label One Side Coated thermoplastic material is a low-temperature thermoplastic material resulting from copolymerization based on polycaprolactone and polyurethane, coated on one side with a foam layer of 0.6 mm thickness. This unique foam layer creates an air layer between the patient's skin and the thermoplastic.

The White Label One Side Coated thermoplastic material hardens within minutes, and adjustments can be made easily. Each sheet is individually packaged and available in various sizes, thicknesses, and types of perforations.

Available in the following sizes:

Article number	Type	Color	Dimensions in cm.	Thickness in mm.
WL/0800	Non-Perfo	Skin color	60 x 41	0,8
WL/0810	Micro-Perfo	Skin color	60 x 41	0,8
WL/1600	Non-Perfo	Skin color	60 x 43	1,6
WL/1610	Micro-Perfo	Skin color	60 x 43	1,6
WL/2000	Non-Perfo	Skin color	60 x 44	2
WL/2010	Micro-Perfo	Skin color	60 x 44	2
WL/2020	Mini-Perfo	Skin color	60 x 44	2
WL/2030	Multi-Perfo	Skin color	60 x 44	2
WL/2500	Non-Perfo	Skin color	60 x 45	2,5
WL/2520	Mini-Perfo	Skin color	60 x 45	2,5

White Label Comfort

One Side Coated

Coated on one side with a foam layer



Ergonomic design

Thoughtfully designed, our products combine comfort with effective functionality.



Sustainable

Optimal consideration has been given to both people and the environment.



Comfortable

The quality of a product is also reflected in the comfort it provides for the patient.



Available in the following sizes:

Article number	Type	Color	Dimensions in cm.	Thickness in mm.
WL/2530	Multi-Perfo	Skin color	60 x 45	2,5
WL/2500/9	Non-Perfo	Skin color	60 x 90	2,5
WL/2520/9	Mini-Perfo	Skin color	60 x 90	2,5
WL/2530/9	Multi-Perfo	Skin color	60 x 90	2,5
WL/3000	Non-Perfo	Skin color	60 x 45	3
WL/3020	Mini-Perfo	Skin color	60 x 45	3
WL/3030	Multi-Perfo	Skin color	60 x 45	3
WL/3000/9	Non-Perfo	Skin color	60 x 90	3
WL/3020/9	Mini-Perfo	Skin color	60 x 90	3
WL/3030/9	Multi-Perfo	Skin color	60 x 90	3
WL/4000	Non-Perfo	Skin color	60 x 45	4
WL/4020	Mini-Perfo	Skin color	60 x 45	4
WL/4000/9	Non-Perfo	Skin color	60 x 90	4
WL/4020/9	Mini-Perfo	Skin color	60 x 90	4



White Label Ortho



Ergonomic design

Thoughtfully designed, our products combine comfort with effective functionality.



Sustainable

Optimal consideration has been given to both people and the environment.



Comfortable

The quality of a product is also reflected in the comfort it provides for the patient.



The White Label Ortho thermoplastic material is a low-temperature thermoplastic material resulting from copolymerization based on polycaprolactone and polyurethane.

The White Label Ortho thermoplastic material hardens within minutes, and adjustments can be made easily. Each sheet is individually packaged and available in various sizes, thicknesses, and types of perforations.

Available in the following sizes:

Article number	Type	Color	Dimensions in cm.	Thickness in mm.
WLOS160	Non-Perfo	Skin color	60 x 43	1,6
WLOS161	Micro-Perfo	Skin color	60 x 43	1,6
WLOS162	Mini-Perfo	Skin color	60 x 43	1,6
WLOS163	Multi-Perfo	Skin color	60 x 43	1,6
WLOS200	Non-Perfo	Skin color	60 x 44	2
WLOS201	Micro-Perfo	Skin color	60 x 44	2
WLOS202	Mini-Perfo	Skin color	60 x 44	2
WLOS203	Multi-Perfo	Skin color	60 x 44	2
WLOS250	Non-Perfo	Skin color	60 x 45	2,5
WLOS252	Mini-Perfo	Skin color	60 x 45	2,5
WLOS253	Multi-Perfo	Skin color	60 x 45	2,5
WLOS320	Non-Perfo	Skin color	60 x 45	3,2
WLOS322	Mini-Perfo	Skin color	60 x 45	3,2
WLOS323	Multi-Perfo	Skin color	60 x 45	3,2

White Label Ortho



Ergonomic design

Thoughtfully designed, our products combine comfort with effective functionality.



Sustainable

Optimal consideration has been given to both people and the environment.



Comfortable

The quality of a product is also reflected in the comfort it provides for the patient.



Available in the following sizes:

Article number	Type	Color	Dimensions in cm.	Thickness in mm.
WLOS320/9	Non-Perfo	Skin color	60 x 90	3,2
WLOS322/9	Mini-Perfo	Skin color	60 x 90	3,2
WLOS323/9	Multi-Perfo	Skin color	60 x 90	3,2
WLOS400	Non-Perfo	Skin color	60 x 45	4
WLOS402	Mini-Perfo	Skin color	60 x 45	4
WLOS400/9	Non-Perfo	Skin color	60 x 90	4
WLOS402/9	Mini-Perfo	Skin color	60 x 90	4

White Label Pré-cuts



Ergonomic design

Thoughtfully designed, our products combine comfort with effective functionality.




Sustainable

Optimal consideration has been given to both people and the environment.



Comfortable

The quality of a product is also reflected in the comfort it provides for the patient.



Duimspalk

The White Label One Side Coated thermoplastic material is also available in various pre-cut models.

Available in the following sizes:

Article number	Type	Color	Thickness in mm.	Size
WL32S1610	Micro-Perfo	Skin color	1,6	Small
WLO32M201	Micro-Perfo	Skin color	2,0	Medium